

# TRAKnet™ DPM Template Creation Quick Notes

Due to popular demand this is a short list of the design and creation tricks for creating your templates.

## Answer Field Creation

Making a choice field (multi-select): [Value 1,Value 2,Value 3]

- o Separate values with commas
- o Leave no extra spaces after commas
- o Other will always be a default value in the encounter

Making a option field (single-select): [OPTION=Value 1,Value 2,Value 3]

- o Option must be in all caps
- o Leave no extra space after the commas
- o Other will always be a default value in the encounter

Defaulting an OPTION or Choice field: [Value 1,Value 2\*,Value 3]

- o To default one or more than one value, put a star after the value, before the comma
- o With an OPTION field you can only default one value

Creating a Question in a Question: [Value 1,Value 2 ??,Value 3]

- o The ?? denotes an addition piece of info needed for a value and you will be prompted to type it in.
  - Example: [Lesion 1 ?? mm long x ?? mm wide x ?? mm deep, Lesion 2 ?? mm long x ?? mm wide x ?? mm deep, Lesion 3 ?? mm long x ?? mm wide x ?? mm deep]

## Data Field Quick Insert

To quickly enter a patient data field into the template type "patient." (minus the quotes)

- o Make sure you type a period immediately after patient
- o A list of patient data fields will appear as a menu
- o Choose the one you would like and the data field will be created in correct format
  - Example: [patient.name]

To quickly enter a provider data field into the template type "provider." (minus the quotes)

- o Make sure you type a period immediately after provider
- o A list of provider data fields will appear as a menu
- o Choose the one you would like and the data field will be created in correct format
  - Example: [provider.name]

To quickly enter a user data field into the template type "user." (minus the quotes)

- o Make sure you type a period immediately after user
- o A list of user data fields will appear as a menu
- o Choose the one you would like and the data field will be created in correct format
  - Example: [user.name]

To quickly enter a referring provider data field into the template type "referral." (minus the quotes)

- o Make sure you type a period immediately after referral
- o A list of referring provider data fields will appear as a menu
- o Choose the one you would like and the data field will be created in correct format
  - Example: [referral.name]